

## Aidacare FL250 Floorline Bed



- Under Bed Clearance - While the FL250 achieves a low height of 110 mm, full clearance under the bed is still maintained, minimising the risk of foot or cord impingement
- SmartDescent™ Automatically slows descent speed over the final 300 mm of lowering to the floor level, allowing possible impingements to be identified and moved
- Hidden Bed Mechanisms: The bed mechanisms are attached under the platform and are not visible. This helps the room feel less clinical
- Low Height 11 cm: Adding a 15 cm mattress gives you a total height below 30 cm thus significantly reducing out-of-bed fall injury risks
- Auto-Regression on back & knee rests: Reduces skin shear forces during operation thus improving clinical outcomes and reducing the need for expensive mattress solutions.
- High SWL 250 kg & base-extension to 100 & 200 mm: The same bed platform can meet a wider profile of residents
- BrakeFirst™ Automatic Safety Brake Activation at lowest height: Reduces the risk of injury during bed exit in the morning as the bed will lock when at its lowest height
- Single Button Seat Position: Easy transfer of weight from lower back to bottom for pressure-care and seated positioning also helps improve lung function for resident
- ExitSafe™ Single Button Ingress/Exit height: Your staff don't have to think about what height to set the bed. Via one button the vast majority of needs are met
- Attendant Control Panel - Allows the caregiver to lock out certain functions and also ensures bed operation is not just restricted to the handset
- Grey Colour Scheme & Narrow-profile frame: Keeps the focus on bedding and furnishings/finishings rather than the bed-frame

✔ Contracted DVA product

✔ Contract approved

### Compliance

IEC 60601-1:2005

IEC 60601-2-52:2005

### Accessories

Click on the cart icon to display the range of optional accessories

- Underbed light standard: Helps to reduce incidence of resident and caregiver falls
- Rubber & Nylon castors: Castors remain on the ground during braking, thus not damaging floors
- King Single Extension Kit (optional retrofit): Allows you to change the profile of your beds (with minimal cost) to meet resident profile changes
- KTL coated anti-corrosion steel: Improved durability for lower lifecycle cost
- Low-maintenance Braking System: Fewer headaches for your MO's!

## Warranty Information

- 10 Years Frame
- 3 Years Electrics
- 1 Year Handset

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
BEB046100	Bed - Aidacare FL250 Floorline - Single - 4 Function - 250kg (Bed ends sep) - Incl: UBL & Extension	2000 mm	900 mm	110 - 740mm	250 kg	105 kg
BEB046150	Bed - Aidacare FL250 Floorline - King Single - 4 Function - 250kg (Bed ends sep) - Incl: UBL & Extension	2000 mm	1070 mm	110 - 740mm	250 kg	115 kg