

Ceiling Hoist - Liko - Multirall 200



- The Multirall 200 hoist is a portable ceiling hoist that can also be used for room-to-room patient transfers
- Can be used to transfer a patient between rooms in a sitting position, easily moving from one rail system to another with the use of a transfer strap
- Its portability makes it an excellent and economical solution for nursing home and hospital care situations where lifting needs are temporary
- Can be easily attached or removed with no tools and minimum effort from ceiling rails to facilitate convenient transport between rooms
- Quick sling bar change thanks to Q-link and Quick-release Hook
- Can be mounted two ways for ease of use

✔ Contracted DVA product

✔ Contract approved

Compliance

AS/NZS ISO 10535:2011

Warranty Information

- Frame: 10 Years
- Electrical Parts: 2 Years

| Code | Display Name | Total Height | Safe Working Load |
|-----------|----------------------------------------|---------------|-------------------|
| LIC359100 | Multirall 200 - Portable Ceiling Hoist | Up to 1600 mm | 200 kg |